

Using Medical Objects: Restarting Word Plugin



Medical Objects Word Plugin will run as a Tray icon and should be visible in the System tray located at the right hand bottom of the desktop.




1. If the Medical Objects Word Plugin tray icon is not visible in the System tray, left click on **"Start"**, select **"All Programs"**. Open **"Start Up"** from the programs listed and select **"Medical-Objects Word Plugin"**.




2. After approximately 20 seconds the Medical Objects Word Plugin icon will reappear in the System Tray.



3. Right click on the Medical Objects tray icon and select **"Login"**.



4. Type the Medical Objects Word Plugin passphrase and left click **OK**. **Note:** The Medical Objects Word Plugin passphrase is usually created at the time of installation. To change the Medical Objects Word Plugin passphrase please see ***Change Medical Objects Word Plugin Passphrase*** guide.



5. Medical Objects Word Plugin is now ready for use.